



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ZNF226-zinc finger protein 226 Gene

Reference: PTXBC017786

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

zinc finger protein 226; Kruppel-associated box protein

For research use only.

GENE DESCRIPTION	zinc finger protein 226
GENE ID	7769
ORIGIN SPECIES	Human
ACCESSIONS	BC017786
VECTOR	pDONR223
NCBI SYMBOL	ZNF226
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatagtgtcaaggaagcagtgacctcaaggacgtggctgtggccttcacggaggaggaattggggctgctgggacctgccagaggaaagctgtaccgagatgtgatgggtggagaaccttaggaacctgctgtagtggggcatccaccttcaacaagatgatcacatagaaa
gaaatgacagcttttgataatgacgacagacccgaagacagggaattatgataccttactgtaaaagctctttgctctatgacctggcctcaaaccttaa