

ARPP-21-cyclic AMP-regulated phosphoprotein, 21 kD Gene

Reference: PTXBC017805

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARPP-21; R3HDM3; RCS; TARPP; cAMP-regulated phosphoprotein 21; R3H domain containing 3; cAMP regulated phosphoprotein 21kDa; cyclic AMP-regulated phosphoprotein, 21 kD; thymocyte cAMP-regulated phosphoprotein; cAMP regulated phosphoprotein 21		

For research use only.

GENE DESCRIPTION	cyclic AMP-regulated phosphoprotein, 21 kD
GENE ID	10777
ORIGIN SPECIES	Human
ACCESSIONS	BC017805
VECTOR	pDONR223
NCBI SYMBOL	ARPP-21
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgagcaaggagacctgaatcaggcaatagcagagggaaggaggactgagcaggagacggccactccagagaacggcattgttaaatcagaaagtctggatgaaggaggagaaactggaactgcagagcggtggaggtcagaatcaagaaagaagaaatccaagtcag
gagcagagaaaaggtaaaactgactcgactctgtgtctgtgaggaaatcttgcacaccaggaggtgaaagtcttcaggatcagactctctga