

FGL2-fibrinogen-like 2 Gene

Reference: PTXBC017813

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	T49; pT49; fibrinogen-like protein 2; fibrinogen like 2		
For research use only.			

GENE DESCRIPTION	fibrinogen-like 2
GENE ID	10875
ORIGIN SPECIES	Human
ACCESSIONS	BC017813
VECTOR	pDONR223
NCBI SYMBOL	FGL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagctggcctaactggctgagctcagctgttttgcgcgttacgggtttgtgtggcaacaatgaacagaggaaattaaagatgaaagagcaaggatgtctgccagtgagactagaagcagagggaatgcgaagaggcaggggagtgccctaccaggttaagc
ctgcccccttgactattcagctcccgaagcaattcagcaggatcaggaggtgttcaagaagtcctcaaaacctcaaggaaatcgtaaatagtctaaagaaatcttccaagactgcaagctcagggctgatgacaacggagaccaggcagaaacggactgtgttaccagtagcaga
gccccgggagaggtgtgtgataacagagttagagaattagagagtgagggttaacaagctgtcctctgagctaaattcatcttctga