



Reference: PTXBC017842

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAP; Kdap; NAP1; NAPA; SNAPA; napsin-A; ASP4; TA01/TA02; asp 4; aspartyl protease 4; kidney-derived aspartic protease-like protein; napsin-1; pronapsin A; napsin A aspartic peptidase		

For research use only.

GENE DESCRIPTION	napsin A aspartic peptidase
GENE ID	9476
ORIGIN SPECIES	Human
ACCESSIONS	BC017842
VECTOR	pDONR223
NCBI SYMBOL	NAPSA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr