

CCL11-chemokine (C-C motif) ligand 11 Gene

Reference: PTXBC017850

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SCYA11; chemokine (C-C motif) ligand 11; eosinophil chemotactic protein; eotaxin-1; small inducible cytokine subfamily A (Cys-Cys), member 11 (eotaxin); C-C motif chemokine ligand 11		

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) ligand 11
GENE ID	6356
ORIGIN SPECIES	Human
ACCESSIONS	BC017850
VECTOR	pDONR223
NCBI SYMBOL	CCL11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggtctccgagcactctgtggtctgtctcatagcagctgcttcagccccaggggtcgtggccagctctgtcccaaccactgctgtttaacctggccaatagggaagatacccttcagcgactagagagctacaggagaatcaccagtggcaaatgtccccagaaagctgtgatcttcaagaccaaactggccaaggatctgtgcccacccaagaagaagtggtgagggattccatgaagtatctggacaaaatctccaactccaaagccataa