

GYPE-glycophorin E Gene

Reference: PTXBC017864

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GPE; MNS; MiIX; glycophorin-E; glycophorin E (MNS blood group)		
For research use only.			

GENE DESCRIPTION	glycophorin E
GENE ID	2996
ORIGIN SPECIES	Human
ACCESSIONS	BC017864
VECTOR	pDONR223
NCBI SYMBOL	GYPE
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtatggaaaaataatctttgtattactattgtcaggaattgtgagcatatcagcatcaagtaccactgggtgtggcaatgcacacttcaacctcttcttcagtcacaaagagttacatctcatcacagacaaatgggataacactcattaattgggtggcgatggctcgtgttattttgaggtgatgct
tgttgtgtgtggaatgatcatcttaatttcttactgtattcagatga