

**FCGR3A-Fc fragment of IgG, low affinity IIIa, receptor (CD16a) Gene**

Reference: PTXBC017865

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD16A; FCG3; FCGR3; FCGR3III; FCR-10; FCRIII; FCRIIIA; IGFR3; IMD20; low affinity immunoglobulin gamma Fc region receptor III-A; CD16a antigen; Fc fragment of IgG low affinity IIIa receptor; Fc fragment of IgG, low affinity III, receptor for (CD16); Fc fra		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	Fc fragment of IgG, low affinity IIIa, receptor (CD16a)
<b>GENE ID</b>	2214
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC017865
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	FCGR3A
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

atgtggcagctgctcctccaactgctctgctacttctagtttcagctggcatgaggactgaagatctcccaaggctgtggttctctggagcctcaatggtacagggctctgagaaggacagctgtgactctgaagtgccagggagcctactcccctgaggacaattccacacagtggttc  
acaalagagcctcatctcaagccaggcctcgagctacttcaatgacgctccacagtcgacgacagtgagagtagcaggtgccagacaacctctccaccctcagtgacccgggtgacagctagaagtcctatcggctggctgctccagccctcgggtgttcaaggagggaag  
accctattcactgaggtgcacagctggaagaacactgctctgataaggtcacatattacagaatggcaaggcagggaagtatttcatcataattctgactctacattccaaaagccacactcaagacagcggctcctactctgcaggggctgttgggagtaaaatgtctctcag  
agactgtaacatcacatcactcaaggtttgagcagtgcaacatctcatcattctccacctgggtcccaagctcttctgctgtggtgatgctactcttttgcagtgacacagactatattctctgtgaagcaaacattcgaagctcaacaagagactggaaggaccataaattaat  
ggagaaggaccctcaagacaatga