

RPP40-ribonuclease P/MRP 40kDa subunit Gene

Reference: PTXBC017871

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RNASEP1; bA428J1.3; ribonuclease P protein subunit p40; RNase P subunit 1; RNaseP protein p40; ribonuclease P (40 kD); ribonuclease P, 40kD subunit; ribonuclease P/MRP 40kDa subunit; ribonuclease P1; ribonuclease P/MRP subunit p40		

For research use only.

GENE DESCRIPTION	ribonuclease P/MRP 40kDa subunit
GENE ID	10799
ORIGIN SPECIES	Human
ACCESSIONS	BC017871
VECTOR	pDONR223
NCBI SYMBOL	RPP40
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacactggaccctattactttgtgaagaattacctcttcataaattacacctgaattcatcagctacattataaagaagggaattattgtcactggataaagacacttatgaagaactggactcagggtcatctcagttttctggcagaaaaattatgaaattattgttccat
tgatttgatggaattatccttaaaccttgattctaaagaagtatgaaagaattcttggtcttcaagaagaagccattgaaatttgattttctttggcctggcataaaacagggtcagaagaatcgacaatgatgcatattttccaagtacaaattcaggagcatcagccaaaagtagcactg
agcacgttgagagatctccagtgccagtgctgcagagcagcgagctggagggaacgccagaggtgtctctgccggctctggagctctcagctggcggcgtcttcagtaattgacgacctaataatgagcctaataattcatacaacctattgctgtcctgagccaagcacagt
ggtggcaaaagcttattgtgtacaatcactggcttcatactccagagaagatgtctctctattggaacatctctgtcactattgatgaaccgaagttagctccatgggttacactgtccgttcaaggctttgcagacagccctgtttctgggaaaaaaatga