

FAM33A-family with sequence similarity 33, member A Gene

Reference: PTXBC017873

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FAM33A; spindle and kinetochore-associated protein 2; family with sequence similarity 33, member A; spindle and KT (kinetochore) associated 2; spindle and kinetochore associated complex subunit 2		

For research use only.

GENE DESCRIPTION	family with sequence similarity 33, member A
GENE ID	348235
ORIGIN SPECIES	Human
ACCESSIONS	BC017873
VECTOR	pDONR223
NCBI SYMBOL	FAM33A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcgagggtcgataagctggaactgatgtccagaaagctgagctgatctggattacattcaatacagctggaatatgaaatcaagactaatcatctgattcagcagtgagaaaaatccagttacactcttaagggaattgtcagtgataaagctcgatatcaaaactttgtatgccgctttaaacagttgctgttgagcagaaagagagtaagagccgatttggctactgtgaaaaagactatgaatatgatacaaaaaactacagaagcaaacagacctggagctgtcaccactgactaaagaagagaaaaactgcggcagagcaattcaaatccatgccagatttatga