

MRPL46-mitochondrial ribosomal protein L46 Gene

Reference: PTXBC017883

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LIECG2; P2ECSL; 39S ribosomal protein L46, mitochondrial; L46mt; MRP-L46; mitochondrial ribosomal protein L46		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L46
GENE ID	26589
ORIGIN SPECIES	Human
ACCESSIONS	BC017883
VECTOR	pDONR223
NCBI SYMBOL	MRPL46
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgcccgtaaggcggacgctgttaggggtggcgggggttggcggcggttcgagagctctggccggcagctctaaagctctgcagcctgttgcagcggcacctcaagcaacggatcccatggcgcttggggcggttgcctgcagcggccacctgtagtct
ccaagccgttgacccattgcaggaagagatggcgtctctactgcagcagattgagatagagagaagcctgtattcagaccacgagcttcgtctctggatgaaaccagcgactggcaagaagaaagctgacctcatgatgaagaagatgaacaggatattgtggcgcaagattt
ggaagatatgtgggagcagaattctacagttcaacttggagctgcataacagaagctgatgaaagaatgaccgaacatccctgaacaggaagctagacaggaacctgtcctgttagtcagagagaagtttggagaccaggatgttggatactgccccaggcagagtggcagcct
ggggagacccttcgaggaacagctgaacgaactctggccacactctcagaaacaacatggaagccaagtcttaggaaatgcaccctgtgggcactacacattcaagttccccaggcaatgcggacagagagtaacctcggagccaaggtgttcttctcaaaagcactgtctattaactg
gagacttttccaggctgggaataaggcccatcatgtgtgggtcactaaggatgagctgggtgactatttgaaccaaataacctggccaagttaggaggtttgttcagacctctga