

SNRNP27-small nuclear ribonucleoprotein 27kDa (U4/U6.U5) Gene

Reference: PTXBC017890

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	27K; RY1; U4/U6.U5 small nuclear ribonucleoprotein 27 kDa protein; U4/U6.U5 tri-snRNP-associated 27 kDa protein; U4/U6.U5 tri-snRNP-associated protein 3; U4/U6.U5-27K; small nuclear ribonucleoprotein 27kDa (U4/U6.U5); small nuclear ribonucleoprotein, U4/U		

For research use only.

GENE DESCRIPTION	small nuclear ribonucleoprotein 27kDa (U4/U6.U5)
GENE ID	11017
ORIGIN SPECIES	Human
ACCESSIONS	BC017890
VECTOR	pDONR223
NCBI SYMBOL	SNRNP27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtcgcagtcgcagccgctctccacggagggaacgtaggcgttcccgggtccacatcccgggagagagaacgcagggcggcagaaaaggtccaggtctcgggagagagatcggagaaggagccgctcgcgatccccgcaccgaagacgtcccgatctccaagacacatag
atccacatctcttccccttctcgaactgaaagaaagaagagatgaggaaagaaagaaacaaagagcaaaagacggcagattactgaggagacttagaggcgaacacagaggagaaatagaaatgatgaagttaattgggatttgcctctttgactccacaaaaggta
agaaggtggatggctctgtaaatgcctatgccataatgtctctcagaagagggaagtacaggcagtcacatgaatcgaaaaggtggattcaacagaccttggatttcattgcatga