

C14orf156-chromosome 14 open reading frame 156 Gene

Reference: PTXBC017895

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf156; DC50; PD04872; SRA stem-loop-interacting RNA-binding protein, mitochondrial; SRA stem-loop interacting RNA binding protein		

For research use only.

GENE DESCRIPTION	chromosome 14 open reading frame 156
GENE ID	81892
ORIGIN SPECIES	Human
ACCESSIONS	BC017895
VECTOR	pDONR223
NCBI SYMBOL	C14orf156
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccctcagcgcgcgaggtgctgcggcgtgcgtagaagtatcaatcagccggtgtctttgtgagaagaattccttgactcggcgctcagtcagctgaaagaacactttgcacagttcgccatgtcagaaggtgcattttaccttttgacaaggagactggcttcacagaggt
ttgggtgggtcagttttcttcagaagaaggacttcggaatgcactacaacaggaaaatcatattatagatggagtaagggtcagggttcacactagaagcgcaaaactccgcaaacatctgatgatgaaaaaagattttga