

PLSCR1-phospholipid scramblase 1 Gene

Reference: PTXBC017901

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MMTRA1B; PL scramblase 1; ca(2+)-dependent phospholipid scramblase 1; erythrocyte phospholipid scramblase		

For research use only.

GENE DESCRIPTION	phospholipid scramblase 1
GENE ID	5359
ORIGIN SPECIES	Human
ACCESSIONS	BC017901
VECTOR	pDONR223
NCBI SYMBOL	PLSCR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacaacaaactcacagatgaatgcttccacggaaacaaactgccagttgggtatctctcagtagtaccgacagcattccaaggacctccagtagatagtggtaccctgggccccaggtcagctaccacccccaccagccggccattcaggtcctggcccagctggct
ttcctgtcccaatcagccagtggtataatcagccagtagatataatcagccagttggagctgcaggggtaccatggatgccagcgccacagcctccattaaactgtccacctggattagaatatttaagtcaggtatttcaagacacaaaatactataaaaaacagaactgtgcttccagcttgc
ttaaccagattagcaaatga