

PDLIM5-PDZ and LIM domain 5 Gene

Reference: PTXBC017902

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ENH; ENH1; LIM; PDZ and LIM domain protein 5; enigma homolog; enigma-like LIM domain protein; enigma-like PDZ and LIM domains protein; PDZ and LIM domain 5		

For research use only.

GENE DESCRIPTION	PDZ and LIM domain 5
GENE ID	10611
ORIGIN SPECIES	Human
ACCESSIONS	BC017902
VECTOR	pDONR223
NCBI SYMBOL	PDLIM5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcaactacagtgtgtcactgtgtgcccagctccttgggttccggctgcaggcggttaaggattcaacatgcctctgacaatctctagtctaaaagatggcggcaaggcagcccaggcgaatgaagaataggcgatgtgtgtctcagcattgatggaataatgcacaaggaaat
gactcatcttgaagccagaataagattaaagggtgtacaggctcttgaatatgactctgcaagagcatctgctgcacccaagcctgagccggttctgttcaaaagaaaacacaagtgacaaataaccctggcactgtgaaatcccacctaagcccaccaagaaaacacattgtga
gcgtatacagagtttatcatgtaccactcacagtgtgcccagcaagaagagactgattgaggatactgaagactggcgtccaaggactggaacaactcaatctgctcttccgaatccttgcccagatcactgggactgaacatttgaagaatctgaagccgataatacaagaaggc
aaaggaaaagatacccttcacgtctttagtcccaatacacaaaattacgtgactggcaccatgaagtttcagcacgtgctcttaacgtacagtga