

TWIST2-twist homolog 2 (Drosophila) Gene

Reference: PTXBC017907

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AMS; BBR SAY; DERMO1; FFDD3; SETLSS; bHLHa39; twist-related protein 2; class A basic helix-loop-helix protein 39; dermis-expressed protein 1; twist basic helix-loop-helix transcription factor 2; twist homolog 2; twist-related bHLH protein Dermo1; twist fam		

For research use only.

GENE DESCRIPTION	twist homolog 2 (Drosophila)
GENE ID	117581
ORIGIN SPECIES	Human
ACCESSIONS	BC017907
VECTOR	pDONR223
NCBI SYMBOL	TWIST2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggagggtccagctgccccgtgccccgtgacagcctgggcaccagcagaggagctcgagaggcagcccaagcgttcgcccgaagcggcgctacagcaagaagtcgagcgaagatggcagccccgacccgggcaagcggcgaagaaggcagccccag
cgcgagctcctcgaggagctgcagagccagcgcatcctggccaacgtgcgcgagcggcagcgccacagctcgtcaacgaggccttcgcccgtcgcgaagatcatccccacgctgccccctgacaagctgagcaagatccagacgtcaagctggccgcaaggtacatagact
tcctctaccaggtcctgcagagcgacgagatggacaataagatgaccagctgcagctacgtggccccagcggcctcagctacgcctctcctggtggcgcagtgaggggcgctggtccatgtccgctccactag