

C1orf64-chromosome 1 open reading frame 64 Gene

Reference: PTXBC017946

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C1orf64; ERRF; ER-related factor; chromosome 1 open reading frame 64		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 64
GENE ID	149563
ORIGIN SPECIES	Human
ACCESSIONS	BC017946
VECTOR	pDONR223
NCBI SYMBOL	C1orf64
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccccgtcagaatcccaggactggagagccaacctcaaaggcaccatccgtgagacaggctggagaccagctccgggtgggaagctggctggccatcagaagaccgtccccacggctcacctgactttgtattgactgcacccacgggaagcagctctccctggcagca
accgcatcaccaccccaagccccagtcaccaatcagggcgtgtcaccacccaatgaagacttacatcgtgttctgtgggaaaactggcccatctgactcgggtgacccccatgggtggggatgccttgcccaggccagggccaccctgcgcctctgcagagggtctgtgcctc
agcttcttcccagtcagccgctctgccccaggaggtcccgaggctaagggaaccctggaaggctgcgcctgtgaggtcttcaacttgggaacagtcaggactcactgaaagccctctcttctgtctgtgggcaggccgattag