



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C2orf47-chromosome 2 open reading frame 47 Gene

Reference: PTXBC017959

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C2orf47; m-AAA protease-interacting protein 1, mitochondrial; matrix AAA peptidase interacting protein 1		

For research use only.

GENE DESCRIPTION	chromosome 2 open reading frame 47
GENE ID	79568
ORIGIN SPECIES	Human
ACCESSIONS	BC017959
VECTOR	pDONR223
NCBI SYMBOL	C2orf47
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcgtggcgcgtctgttgcaccacagcttctgcactcgcgtcgtcgcctcggggcgcctcgactcgcactcgcctcgtcgtcggcagagtgagagctgcctcggccacacttctgactctcgcctcgtcgcctcgccttgggagcggcggtatttccacgaagcgtatggcctt
ggcagcctcggcgctactcgtccacgggctccgggtggccagtgctcagcagcccgggactccccgagccttcgtcttcttccctcctcgcctcagcgcagctacagcacggaggagaagccccagcagcaccagaaaccaagatgatcgtctcgggattctcaaccccatcaac
tgggttaggactcgaattaaagcctcttacttgcgcctatttgcgaaagagtgatcagatcacagaggttctccagaggagcgaagcagcgtctcgtcatgtatccaagtgtcgtcagagtgtaaaatttgatctgtggaagaacttctgcccgaagaggtgctcatagcatgaagaaaaag
gttacttctactactgaacacataaaatgccctctgcctctcaactaatgatattgaactatgacacagagacatcttacttactatgatgagaagaaggaagtgttgaatacctcgtctgtcttttggatcttaccagcttgcgaacatccccagtcggaactgtga
ttccagtaattggtgggaacacagaatttgaacataaacacattcttatgcaactgtatttgcagaggagttcacacaaagcgtgaagacgtcactgaacatttgatgaactcgaattatagaataa