

HLA-DPB2-major histocompatibility complex, class II, DP beta 2 (pseudogene) Gene

Reference: PTXBC017967

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	major histocompatibility complex, class II, DP beta 2 (pseudogene)
GENE ID	3116
ORIGIN SPECIES	Human
ACCESSIONS	BC017967
VECTOR	pDONR223
NCBI SYMBOL	HLA-DPB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgatcctgcaggtttcagggggccctggacagtggctctgacagcattactgatgggtgctgctcattctgtgtccagagcagggccactccagagaattccgtctaccaggaacggcaggaatgctatgcttcaatgggactcagcgcgttggacgggctcatctacaaccg
ggaggaatacgtgcatttgacagcgcagtgggggagttcctagcagtgatggagctggggcggccatagcgcagtaactcaatagccagaaggactttatggaacggaagcgagccgaggtggacaaggtgtgcagacacaagtacgagctgatggagccactcatccggcagc
gccgaggagacgtaaccataactgctgttagggggtgttggacgacgattctttctgctactctctgtgaaaaggggcgtcgtgtcggggggcgtgcagttggggctcctcctga