

C1orf162-chromosome 1 open reading frame 162 Gene

Reference: PTXBC017973

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein C1orf162; chromosome 1 open reading frame 162		
For research use only.			

GENE DESCRIPTION	chromosome 1 open reading frame 162
GENE ID	128346
ORIGIN SPECIES	Human
ACCESSIONS	BC017973
VECTOR	pDONR223
NCBI SYMBOL	C1orf162
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaggcaatggctccacatgtaaacccgacactgaaagacaaggcactctctccacagcagccccaactagccctgcacccgtgtctcttaaccaccacaacaaaaacattaatccttgccctttgtgtgggttctactgacactgctgctgatagcccttatctctcatcata
aagagctacagaaaatacgaaggagaggttcccaatccccagggcccatgctgagatgggtccctctcttctgacacatggcagatcactccaagcccgagcccgagatcctcactcagcccaagctttcatccatccaggggaatcacttacctatgccagcaca
actttcaactctcagaagaaaagagcaatcacttggctgagaacattctgcagacttggacccattgtctatgctcaattaaagtaacaaactaa