

C8orf22-chromosome 8 open reading frame 22 Gene

Reference: PTXBC017981

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	pancreatic progenitor cell differentiation and proliferation factor-like protein; exocrine differentiation and proliferation factor-like protein; chromosome 8 open reading frame 22		

For research use only.

GENE DESCRIPTION	chromosome 8 open reading frame 22
GENE ID	492307
ORIGIN SPECIES	Human
ACCESSIONS	BC017981
VECTOR	pDONR223
NCBI SYMBOL	C8orf22
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcatcgtaccttccattgggtgccttctagccagaaatcagattatcgaagtcagtggtttctcagttgctttaaactagctctgattctgttaactcatagatgacgacaaaccacagcaaggggtgcctgaagtggcagaatccacctgggtgttaaatggttttcattctgaacctgtgctttcaaatgtgagaataaaagatctgtctgctactggcatgaacagttga