

H2AFZ-H2A histone family, member Z Gene

Reference: PTXBC018002

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2A.Z-1; H2A.z; H2A/z; H2AZ; histone H2A.Z; H2AZ histone; H2A histone family member Z		

For research use only.

GENE DESCRIPTION	H2A histone family, member Z
GENE ID	3015
ORIGIN SPECIES	Human
ACCESSIONS	BC018002
VECTOR	pDONR223
NCBI SYMBOL	H2AFZ
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctggcgttaaggctggaaggactccggaaggccaagacaaaggcgtttccgctcgcagagagccggttcagttccagtgggccgtattcatcgacacataaatctaggacgaccagtcattggacgtgtgggcgcgactgccgctgtgtacagcgcagccatcctgg
agtacctaccgcagaggtactgaaactggcaggaatgcatcaaaagacttaaggtaaaagcgtattaccctcgtcacttgcaacttgctattcgtggagatgaagaattggattctctcatcaaggctacaattgctggtgggtgtcattccacacatccacaatctctgattgggaaga
aaggacaacagaagactgtctaa