

C10orf4-chromosome 10 open reading frame 4 Gene

Reference: PTXBC018007

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C10orf4; F26C11.1-like; FRA10A; protein FRA10AC1; fragile site 10q23.3; fragile site, folic acid type, rare, fra(10)(q23.3) or fra(10)(q24.2) candidate 1; rare folic acid-type fragile site, FRA(10)(q23.3), candidate gene 1; FRA10A associated CGG repeat 1		

For research use only.

GENE DESCRIPTION	chromosome 10 open reading frame 4
GENE ID	118924
ORIGIN SPECIES	Human
ACCESSIONS	BC018007
VECTOR	pDONR223
NCBI SYMBOL	C10orf4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcatggcatggagctatgattctgattttatgtgatgaacactgtggagaatccagcaaaaggaaaaaggacagttgaagacgacttactgtccaaaaaccatttcagaaagaaaaacatggaagggtggctcataaagttgcagcagaattgctgatagggaagaagc
aagaaatagaaggtttcatctcatagctatggatcttatcaagacatagaagttcgttaaatgactatattttatactatgttgcaaaaaaagacttcaagcgtttgggggaaaaatgacaagacagacttggatgtttatagcaaaaaatcatagattcctatggaatgaggaggacgaat
ggacatgacttgggagaagagacttgcataaatactatgataaattatttaagggaatactgcatagcagatctcagtaaatataaagaaaaaagttggatttaggtggcgagtagaaaaagaagtaatttcaggaaaaaggtcaattttctgtggaataaataattgtgataaaaaagagc
ttaaagagttgggaagttaatttggttatattgagcatggtgagaagagaatgcactgttaataaagggttatgccagaatgttccattaaatttaattccatcacaggagaaaaaagaaatcaagtcaaaaaaaagaaaagataaaacaaaaagactgtgaagagtcacataaaaa
tccagattatctctgcagaagaggcctccaagaaaaagataaaggacattctcaagaagaatcgaagattctctacttagaaactctgatagggaagaaagtcctcagaatctgaacttggaaagggtccactaccagagacagatgaaaaatcacaggagaagaatttgatgag
atttcaggattgtttctatga