

FOXP2-forkhead box P2 Gene

Reference: PTXBC018016

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CAGH44; SPCH1; TNRC10; forkhead box protein P2; CAG repeat protein 44; forkhead/winged-helix transcription factor; trinucleotide repeat containing 10; trinucleotide repeat-containing gene 10 protein; forkhead box P2		

For research use only.

GENE DESCRIPTION	forkhead box P2
GENE ID	93986
ORIGIN SPECIES	Human
ACCESSIONS	BC018016
VECTOR	pDONR223
NCBI SYMBOL	FOXP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgcaggaatctgcagacagacaataagcaacagttcaatgaatcaaaatggaatgagcactctaagcagccaattagatctggcagcagagatggaagatcaagtggtagaccagctctgaagtaagcacagtagaactgctcatctgcaacaacagcaggagatgtgtc
tcttacaccagggtattttgtaa