

TSPAN7-tetraspanin 7 Gene

Reference: PTXBC018036

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	A15; CCG-B7; CD231; DXS1692E; MRX58; MXS1; TALLA-1; TM4SF2; TM4SF2b; tetraspanin-7; CD231 antigen; T-cell acute lymphoblastic leukemia associated antigen 1; cell surface glycoprotein A15; membrane component chromosome X surface marker 1; membrane componen		

For research use only.

GENE DESCRIPTION	tetraspanin 7
GENE ID	7102
ORIGIN SPECIES	Human
ACCESSIONS	BC018036
VECTOR	pDONR223
NCBI SYMBOL	TSPAN7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcatcgaggagaatggagacaaacctgtgataacctgtctaaaacctctcatcatctactctctctctctggaactggtgatctgctggtgttgagctctggggcaaaactactctgggacatatctccctattgccaagaactccaaatgctccctatgtgctcat
cggaactggcaccactattgtgtcttggcctgttggatgttgtacatgtctgtgtagccatggatgctgaaactgtatgccatgtttctgtccctgtgttctgtgctgagctgtagctggcatttcagggtttgtgttctgcatgagatcaaggacaccttctgaggactacacggac
actatgcagactacaatggcaatgatgagaggagccggcagtgacatgtgcagcgagcctgagctgtgtgtgtgcagaactacaccaactggagcaccagccctacttctggagcatggcatccccccagctgctgcatgaacgaaactgattgtaatcccaggatcta
cacaactgactgtgcccgcaccaaagttaaccagaagggtgttatgatctggtactagtttcagagactaacatgggaatcatgctggagtggcgtttggaatcgacttcccagtaattggcatgctgctggcctgctgtctgtcccgttcatcacggccaatcagtagatg
gtgtaa