

WDR37-WD repeat domain 37 Gene

Reference: PTXBC018044

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WD repeat-containing protein 37; WD repeat domain 37		
For research use only.			

GENE DESCRIPTION	WD repeat domain 37
GENE ID	22884
ORIGIN SPECIES	Human
ACCESSIONS	BC018044
VECTOR	pDONR223
NCBI SYMBOL	WDR37
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccacagaagcgcaagttgttcgactactcgccaaacaaacagagaagcgcaatccatagcctttctatcgaagaactaacagctcggagcaggagagggagcgggactgccaagagacatgttagaaggacaagattctaaactgccttcctcggttcgcagtacactcttgaa
ctgtttgtcaaatagaaagagaatttgaaaacctttatatcgaaaacttagaattacgtagagaaatcgacactttaataacgttttagctgctgaaggacaagcgattgatggagcagagctgagtaaggccaaactcaaaacaaaagccagtcacagcaccagccagctctccagaaa
ctgaagaccacttacaaggcttcaccagcaagattgtctccagctttaagaccacgacatcgagagctgcctgccagctcgtgaaggagatcatcggccaccgggacggcatctgggatgtcagcgtggccaagacacagccagtggtgctcgggactgcatcagccgatcacacgg
ctttgctgtggagcatagagacagggaagtgcttagtcaagtagcagggccacgtgggctcagtaaatctatcaaaatttcatccctcagagcagttggctctcactgcttctggagataagactgctcatatctggagatacgcgggtgcagctgccgacccccagcctgttgacactag
tgtaagcacccttcctacctgcgaatgtgtaggatgctgatgtctcggaatctgagaatccattatga