

PFN2-profilin 2 Gene

Reference: PTXBC018049

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	D3S1319E; PFL; profilin-2; profilin II; profilin 2		
For research use only.			

GENE DESCRIPTION	profilin 2
GENE ID	5217
ORIGIN SPECIES	Human
ACCESSIONS	BC018049
VECTOR	pDONR223
NCBI SYMBOL	PFN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccggttggcagagctacgtggataacctgatgtcgatggctgctgccaggaggccgcatgtcggtactgcgacgccaaatagctctgggcagccacggcagggggcgtctttcagagcattacgccaatagaaatagatgatttaggaaaaagaccgggaaggtttcttta
ccaacggtttggctcttggcgcgaagaaatgctcagtgatcagagatagctatagctcgatggtgactgcacaatggacatccggacaagagtcagggtggggagccaacatacaatgtggctgtcggtagagctggtagagcttggctcttgaatgggaaaaaggggtccatgg
aggcggattgaataagaaggcactactcaatggcacaatacttgagagactctgggttctag