

LOC646470-similar to proteasome (prosome, macropain) 28 subunit, 3 Gene

Reference: PTXBC018064

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	similar to proteasome (prosome, macropain) 28 subunit, 3
GENE ID	646470
ORIGIN SPECIES	Human
ACCESSIONS	BC018064
VECTOR	pDONR223
NCBI SYMBOL	LOC646470
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcttttctgctgaagctagagccggacatgcagggaagtagattgtttcagagccgtatcggcaggagccgaggatctgtgtccacettcttctcagaaggtcttgagctgaatagccgagttcaggagctccgctgcaagacctgtccagaatccattcggtgccaacc
tcggagccccctggccaccaggaacaaggagatgggccaacaaaatctgcagttctgtgactcagtttccatcaaggtccagcactgttgcagcgaggagacagcttctgcggagcaaccagcatctggtggagttgactgagtggtgaaacctgagatcaagctgct
gagagagaaatgcaacacagtgacatgtag