

C11orf73-chromosome 11 open reading frame 73 Gene

Reference: PTXBC018080

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C11orf73; HLD13; HSPC138; HSPC179; L7RN6; OPI10; protein Hikeshi; lethal, Chr 7, Rinchik 6; Hikeshi, heat shock protein nuclear import factor		

For research use only.

GENE DESCRIPTION	chromosome 11 open reading frame 73
GENE ID	51501
ORIGIN SPECIES	Human
ACCESSIONS	BC018080
VECTOR	pDONR223
NCBI SYMBOL	C11orf73
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtttgctgctgtgtggcggggaggctggtgcaacagctgcacagcaagtgccagagataaattgttttgactacctgattatgaaagtatcaacatgttgtgttttatgctgggaacaatccattgctgagggaatgggaggatctgtctactttcttatctgattcaaatgga
atgccagtatggcaactcctaggattgtcacgaatgggaagccaaagtgccatctcaaaatttcaggtctaaatctggagaaggaagccaacatccttttgagccatgaatattgtccgaactccatctgtgctcagattggaatttcagtggaaattatgacagtagctcagcagactc
ctgtaggtaaatgctgtgtatcctcagttgactcattcactcagttcacacaaaagatgttgacaattctacaattttgcttcatcattgtctgtctcaggccagatgacaccaagcccatctgaaatgttcattccggcaaatgtgttctgaaatgtgtatgaaaactttcaagacgactagca
cagaacccctctcttttgaaaacataa