

C6orf57-chromosome 6 open reading frame 57 Gene

Reference: PTXBC018085

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0369 protein C6orf57; C6orf57; Sdh8; succinate dehydrogenase assembly factor 4, mitochondrial; SDH assembly factor 4; succinate dehydrogenase complex assembly factor 4		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 57
GENE ID	135154
ORIGIN SPECIES	Human
ACCESSIONS	BC018085
VECTOR	pDONR223
NCBI SYMBOL	C6orf57
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacccatcgaggctccctgggtgcttagctgggtctcgccacggcgtagagcggcaagatcacccctctgtgctattctctgaggaaaacaagttcttctcaaggaggaaagtctgaactgtcaaacagtccttaagaagccgatgttaccagaaggtcgtttgatgcaccagaggattccatttagagaagaaccactggaaaaattccagatgatgtaatccagtgaaccaagaaaaaggtggaccagggggccagaaacctaccgatatggagattgggaacgaaaaggacgctgtattgatttttaa