

RNF141-ring finger protein 141 Gene

Reference: PTXBC018104

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZFP26; ZNF230; RING finger protein 141; C3HC4-like zinc finger protein; zinc finger protein 230		

For research use only.

GENE DESCRIPTION	ring finger protein 141
GENE ID	50862
ORIGIN SPECIES	Human
ACCESSIONS	BC018104
VECTOR	pDONR223
NCBI SYMBOL	RNF141
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggacagcaaatcggatcagacacagttggttattaacaagtaccagaaaagtagcaaaacatgttacgttgggtcgagagagtggtccttaacttatgaagaatttctgggagagtagctgagcttaatgatgaacggcctaagtggtcttgccaggaaaaacatcttctttg
aggtacaacctgggtctgattcctctgcttttggaaagtgggtgacgggtggtctgtaccaagattaacaaaagcagtggtcatttggaggcatcacggatcatgaatttataccagttattcaactttataaagatacacaaagtaagcagcaggagattggcacagagctccacctctga
agaacctgatgaaaactatcctctgttaacatctgtcaggctagcttttggatgggaagggtgaagcagctgaccgatgaggaggagtggtgtatctgtatggaatggcgggctgacctcatcgccttgctcacagctttgtcagaagtgattgataaatggagtgatcgacacagga
attgccctatttgcctacagatgactggagcaaatgaatcttgggtggtatcagatgcacccactgaagatgataggtaactatattcttaacatggctgatgaggcaggccagccacagggccatga