

RNF138-ring finger protein 138 Gene

Reference: PTXBC018107

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase RNF138; HSD-4; NARF; STRIN; hNARF; NLK-associated RING finger protein; Nemo-like kinase-associated RING finger protein; nemo-like kinase associated ring finger protein; ring finger protein 138, E3 ubiquitin protein ligase; test		

For research use only.

GENE DESCRIPTION	ring finger protein 138
GENE ID	51444
ORIGIN SPECIES	Human
ACCESSIONS	BC018107
VECTOR	pDONR223
NCBI SYMBOL	RNF138
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgaggacacctctcggccacgtctacaccgaagatgattctactgcccgtctgcaggaggtgtctcaaaacgcccgtgcggaccacggcctgtcagcacgtttctgtagaaaatgttctgactgcaatgagggaagcggagcacattgtcccctatgtcgtggaatgtg
actagaagagagagagcatgtcctgaacggccttagaccttgaaaataatgagggaagtttctggtgctgcagatgctgtgcacaaacagattaaattctatcgcatgagacatcattacaatcttgtaagaagtacaggatgaatatgggtgtttctctatcattccaaacttcagatctct
caagattcagtagggaacagcaataggagtgaaacatccacatctgataacacagaaacttaccaagagaatacaagttcttctggtcatcctacttttaagtgtccctgtgtcaagaatcaaatttaccagacagcgtttactggatcactgtaacagtaacacatttcagatagttcctgtg
acatgtcctattgtgtctctccttgggagatcctagccagattaccagaaattcgttagtcatcctaatacagagacatcaattgattatggagaattgtgaattctcagctagatgaagaacccaataccaactcgtgtgaagaattttcaagtaaacatctga