

TPD52-tumor protein D52 Gene

Reference: PTXBC018117

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	D52; N8L; PC-1; PrLZ; hD52; tumor protein D52; prostate and colon associated protein; prostate leucine zipper; protein N8		

For research use only.

GENE DESCRIPTION	tumor protein D52
GENE ID	7163
ORIGIN SPECIES	Human
ACCESSIONS	BC018117
VECTOR	pDONR223
NCBI SYMBOL	TPD52
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaccgcggcgagcaaggtctgtgagaacagaccagtcctgaggaaggagaagatgttctgccacgacagtcacagagaccctctcgaaggaggagcaggagagctaaagaagaactgcaaaaggtagaagaagaatccagactctgtctcaagtgttagcagc
aaaagagaagcatctagcagagatcaagcggaaacttggaatcaattctctacaggaactaaaacagaacatgccaagggtgcaagacgtgacagcaacatctgcttaagaagacatctgaaccttatccaggctggacagaaggcctcagctgtttctgtgttggtcag
tcacacccaaaagctggaagatgtaaaaaactcccaacttttaatactttgaagaaaaggcgaacacaaagtctaaagttagggggaaccaagcctgctgtgtgtgattttggagaagcttgaattcggtgcaaatgtagtgcaccaccacggagcctcttcagaaaaagacac
aggagagcctgtga