

C12orf65-chromosome 12 open reading frame 65 Gene

Reference: PTXBC018145

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COXPD7; SPG55; chromosome 12 open reading frame 65		
For research use only.			

GENE DESCRIPTION	chromosome 12 open reading frame 65
GENE ID	91574
ORIGIN SPECIES	Human
ACCESSIONS	BC018145
VECTOR	pDONR223
NCBI SYMBOL	C12orf65
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcaccgtgggtttatttctacaccactgacccgaatatgcccgccatggggactcggctttgggagaagctgacgtgttatccccaggaatagctgtcactccggtccagatggcaggcaagaaggactacccctgcactgcttcttgatgagaatgaactcgaag
agcagtttgtaaaggacacgggtccaggggcccaggaaccaacaaaaccagcaactgcgtgtgtctgaagcacatccccctcaggcatcgtgttaaagtgccatcagacaagatcagttgatcagaacagaagctagctcggaaaatcctacaagagaagtagatgtttctacaatg
gtgaaaacagtcctgttcacaagaaaaacgagaagcggcgaaagaaaaacaagaaggaaaaaaagagcaaaaggaaacctggaaaaaaagaagctactaaagaactgtgggagtcagtaaaaagggtccactga