

IL15-interleukin 15 Gene

Reference: PTXBC018149

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL-15; interleukin-15; interleukin 15		
For research use only.			

GENE DESCRIPTION	interleukin 15
GENE ID	3600
ORIGIN SPECIES	Human
ACCESSIONS	BC018149
VECTOR	pDONR223
NCBI SYMBOL	IL15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagaattcgaaccacatttgagaagtatttccatccagtgtactgtgtttacttctaaccagtcattttctaactgaagctggcattcatgtcttatttgggtgtttcagtcagggtctcctaaaacagaagccaactgggtgaatgaataagtgatttgaataaattgaagatcttattc
aatctatgcataatgatgtactttatatacggaaagtgtgttcccccagttgcaagtaacagcaatgaagtgctttctctggagttacaagttatttcacttgagtcggagatgcaagtattcatgatacagtagaaaatctgatcatcctagcaaacacagttgtcttctaatgggaatga
acagaatctggatgcaagaatgtgaggaaactggaggaaaaaataaagaattttgcagagttttgtacatattgtccaaatgtcatcaacactcttga