

CFL1-cofilin 1 (non-muscle) Gene

Reference: PTXBC018256

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CFL; HEL-S-15; cofilin; cofilin-1; 18 kDa phosphoprotein; cofilin 1 (non-muscle); cofilin, non-muscle isoform; epididymis secretory protein Li 15; p18; cofilin 1		

For research use only.

GENE DESCRIPTION	cofilin 1 (non-muscle)
GENE ID	1072
ORIGIN SPECIES	Human
ACCESSIONS	BC018256
VECTOR	pDONR223
NCBI SYMBOL	CFL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctccggtgtggtgtctctgatggtgtcatcaaggtgttcaacgacatgaaggtgcgttaagtcctcaagccagaggaggtgaagaagcgaagaaggcgtgtctctctgctgagtgaggacaagaagaacatcatcctggaggagggaaggagatcctggtggcgatgt
gggccaagactgtcgacgacccctacgccactttgtcaagatgctccagataaaggactgccgctatgccctctatgatgcaacctatgagaccaaggagagcaagaaggagatctggtgtttatctctggcccccagctcgccccccttaagagcaaaatgattatgccagctcca
aggacgccatcaagaagaagctgacagggatcaagcatgaattgcaagcaaacctctacgaggaggtcaaggaccgtgcaccctggcagagaagctggggggcagtgccgtcatctcctggagggaagcctttgtga