

## PFKL-phosphofructokinase, liver Gene

Reference: PTXBC018295

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | ATP-PFK; PFK-B; PFK-L; ATP-dependent 6-phosphofructokinase, liver type; 6-phosphofructokinase type B; 6-phosphofructokinase, liver type; liver-type 1-phosphofructokinase; phosphofructo-1-kinase isozyme B; phosphohexokinase; phosphofructokinase, liver type |                     |                |

**For research use only.**

---

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | phosphofructokinase, liver                                                                       |
| <b>GENE ID</b>          | 5211                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC018295                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | PFKL                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggcgatggagcggtgatggcgtgctggaagccacgctgacacgcccgcctgctgtgtcaccctctcggggaaccagtcagtcgggtgccccctcatggagtgcgtgcagatgaccaaggaagtcagaaagccatggatgacaagaggtttgacgaggccaccagctccg  
tgggtggagcttcgagaacaactggaaacatttacaagctctcgcaccagaaagcccccaaggagaagtctaacttccctggccatctgaatgtggggccccggcgtggcatgaatgcggccgtgcgctcggcggtgcggaccgcatcccatggacacagtatac  
gtgtgtcacatggtctcgaagcctagccaaggtcaggtgcaagaagttaggtggcagcagctggccgctgtgtggggcgtgtgtgtccatgctggggaccaagagaccctgcccaggccagctggagtcattgtggagaacatccgcatctatgtattcacgcccgt  
ctgtgtgtcgtgtgtgtgagccctatgaaggggtgctgcagctgtggaggtcgcggcgctacgaggagctctgcatgctatgtgtgtatccagccaccatcagcaaacgtccctggcaccgacttcagcctgggctccgacactgctgtaaatgcccatggagagctgt  
gaccgcatcaaacagctctgctcggggaccaagcgcctgtgtgtcatcgtggagacatgggggttactgtgtgctacctggccaccgtgactggcattgctgtggggccgacgcccgtacgtcttcgaggacccttcaacatccacgacttaagggtcaacgtggagcacatgacg  
gagaagatgaagacagacattcagaggggctgtgtgtcgtgggaacgagaagtgccatgaactactacaccacggagttcctgtacaacctgtactcatcagagggcaaggcgctcttcgactgcaggaccaatgtcctggggccactgcagcaggtggcgctccaacccctttgacc  
ggaaactatgggaccaagctgggggtgaagccatgctgtgtgtgtggaagagctgcgcgaggtttaccgcaaggggacgggtgtgtgcgaatgcccagactcggcctgctgtgctggcctgaagaaagagcggtgtggccttcagccccgtcactgagctcaagaaagacactgattt  
cgagcaccgcatgccacgggagcagtggtgtgtgagcctgctggtcctatgctgaagatgctggcacaataccgcatcagatggcgccctacgtgtcaggggagctggagcacgtgacccgccgacccctgagcatggacaagggtcttga