

NFATC2IP-nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 2 interacting protein Gene

Reference: PTXBC018311

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ESC2; NIP45; RAD60; NFATC2-interacting protein; 45 kDa NF-AT-interacting protein; 45 kDa NFAT-interacting protein; nuclear factor of activated T-cells, cytoplasmic 2-interacting protein; nuclear factor of activated T-cells, cytoplasmic, calcineurin-depend		

For research use only.

GENE DESCRIPTION	nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 2 interacting protein
GENE ID	84901
ORIGIN SPECIES	Human
ACCESSIONS	BC018311
VECTOR	pDONR223
NCBI SYMBOL	NFATC2IP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggagccctgcagagtggtggaccacatggccaccacgtgggtgtcccaagcaggatccttttcttttggagagacagatcacctactgccactccaggaccctaaagctcggagtggctgacatcattgactgtgtgtactaacaagtctccagaggccaca
gagacgtcccaacagctccagctccgggtgcaggaaaggagaaacaccagacactggaagtctcactgtctcagatgccctctaaagacctcatgtcccactatgaggagccatgggactgtcgggacggaagctctctctctttgatgggacaaagcttcaggcaggagg
ctgccagctgacctggcgtggaatctggggacctcattgaggtctggggctga