

RABIF-RAB interacting factor Gene

Reference: PTXBC018488

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MSS4; RASGFR3; RASGRF3; guanine nucleotide exchange factor MSS4; Ras-specific guanine-releasing factor 3; mammalian suppressor of SEC4; rab-interacting factor; RAB interacting factor		

For research use only.

GENE DESCRIPTION	RAB interacting factor
GENE ID	5877
ORIGIN SPECIES	Human
ACCESSIONS	BC018488
VECTOR	pDONR223
NCBI SYMBOL	RABIF
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaccagcggagcagccgagcaggttagtgcagccgagggccgaaaccggaaggcgggtgtgtgcagcgttgcggctccgggtgtgtgcagccagggacgctctctctcgcgcagcgtttcttcctccatgagaagaagccagctctgtgtgacggcagcaatc
ctgacggcgatctctccaggaacactggctggttgaggacatgttcattttgagaatgtgggcttcaccaaggacgtgggcaacatcaagtttctgtgtgcgcagactgtgaaattggaccaattggctggcattgcctagatgacaagaacagtttctatgtggccttgaacgagtttcc
catgagtaa