



Reference: PTXBC018514

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TC1; THMD1; THT1; THTR1; TRMA; thiamine transporter 1; high affinity thiamine transporter; reduced folate carrier protein (RFC) like; solute carrier family 19 (thiamine transporter), member 2; thTr-1; solute carrier family 19 member 2		

For research use only.

GENE DESCRIPTION	solute carrier family 19 (thiamine transporter), member 2
GENE ID	10560
ORIGIN SPECIES	Human
ACCESSIONS	BC018514
VECTOR	pDONR223
NCBI SYMBOL	SLC19A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]