

FAM36A-family with sequence similarity 36, member A Gene

Reference: PTXBC018519

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FAM36A; cytochrome c oxidase protein 20 homolog; COX2 Cox2 chaperone homolog; family with sequence similarity 36, member A; COX20, cytochrome c oxidase assembly factor		

For research use only.

GENE DESCRIPTION	family with sequence similarity 36, member A
GENE ID	116228
ORIGIN SPECIES	Human
ACCESSIONS	BC018519
VECTOR	pDONR223
NCBI SYMBOL	FAM36A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgccccgcggagcccggtgagcccgaggagaggaagtccttaagctcctaggaatttttagatgtgaaatactccctgcgcccgccatcaatattgtatggttcattagatctgtgtggctggctttgacattttttgtcactagtagaattagaagatcatgtgatgttgagtaggaggtttatcttggtagcttgggagctggtttcattgtaggtataattatgcaagcaagaatccaggaagaattgccagagaagaataaaaaagaagataattatgaaggtacccacctgatctgaaagaaaacacaacggcagcagcaattga