

GPR157-G protein-coupled receptor 157 Gene

Reference: PTXBC018691

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	G protein-coupled receptor 157		
For research use only.			

GENE DESCRIPTION	G protein-coupled receptor 157
GENE ID	80045
ORIGIN SPECIES	Human
ACCESSIONS	BC018691
VECTOR	pDONR223
NCBI SYMBOL	GPR157
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagccgtccccgcccccacgagctgggtccgtcggagcggccgtgggtgctgtgtgctgcgcactctccgcgtcgggtcggcctgtgggtggccacgcacgccctgtggcccgacgtgcgcagccgggcacggcgcctgtgtctctctgtgctggccgacctgtct
cggccgctctacttctacggagtgctgcagaacttcggggccgtcgtgggactgcgtgctgcagggcgcgctgtccaccttgcgaacaccagctctctctgtgaccgtggccattgcgtctacttgcacctcagcatgtccgcgccgcgcgggctcgcacagatcgct
gctttggccttccatgtcgtcaggtgggtggcggtggcgctgtttccaggagccccgacacaggccgacccctccggtcttgcctcccagaggcgcgtctag