

LDOC1L-leucine zipper, down-regulated in cancer 1-like Gene

Reference: PTXBC018713

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein LDOC1L; Mar6; Mart6; dJ1033E15.2; leucine zipper protein down-regulated in cancer cells-like; mammalian retrotransposon-derived protein 6; leucine zipper down-regulated in cancer 1 like		

For research use only.

GENE DESCRIPTION	leucine zipper, down-regulated in cancer 1-like
GENE ID	84247
ORIGIN SPECIES	Human
ACCESSIONS	BC018713
VECTOR	pDONR223
NCBI SYMBOL	LDOC1L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgctctttaagcgtgatcgagggaagacacagcaggccaccattccatgaatgggaggtgtacagatcattctcttctgtcagttctgtctgtcctccagcagctatatgtgaagactagtagctgccaggagaggtgcccccaagtgaaggggtacagtggcacctgg
gaaaaggcacctggaaggtttccatgtggccagcccagcatggaagcagggtgggaactctgtgtgtcgcagcgcctcactctactcagtggtcttttgaagccctaccatgtctgtcaggcctgtgctgtcttcacatctacagctgcctag