



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

GPSM3-G-protein signaling modulator 3 (AGS3-like, *C. elegans*) Gene

Reference: PTXBC018724

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

AGS4; C6orf9; G18.1a; G18.1b; G18.2; NG1; G-protein-signaling modulator 3; G-protein signaling modulator 3 (AGS3-like, *C. elegans*); G-protein signalling modulator 3 (AGS3-like, *C. elegans*); activator of G-protein signaling 4; G-protein signaling modulator

For research use only.

GENE DESCRIPTION	G-protein signaling modulator 3 (AGS3-like, <i>C. elegans</i>)
GENE ID	63940
ORIGIN SPECIES	Human
ACCESSIONS	BC018724
VECTOR	pDONR223
NCBI SYMBOL	GPSM3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagagctgagagaccacccaggagaagagatggatggacagggccccctcaggatgagggaaggctggccccctcaaaactccacaactcggcctggcgatctgctctccatccctctctctccaggagaccgccacacagccctgggaccccgctcggcctccctgctctccctcagactgaactctctggactctggctggctggcggccagctcgccgcctggagagagagagggccactctacaccccccaaacccctcaagctccgcctctccgccactctgctctctcaggagagagagaacatctatagacactatctatgacacagctgccagctggatggagccagcctgctgaagctctccctccctccagggcggaagactctctggggcggaagactctctggatctctcagaaatctacagctggggctggaatggagagcaagctgcccgcggcccccaacacactctgtag

More info: Please contact us for any question: contact@proteogenix.fr