

CD74-CD74 molecule, major histocompatibility complex, class II invariant chain Gene

Reference: PTXBC018726

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD74 molecule; CD74 molecule, major histocompatibility complex, class II invariant chain; CD74 antigen (invariant polypeptide of major histocompatibility complex, class II antigen-associated); DHLAG; HLADG; Ia-GAMMA; HLA class II histocompatibility antigen		

For research use only.

GENE DESCRIPTION	CD74 molecule, major histocompatibility complex, class II invariant chain
GENE ID	972
ORIGIN SPECIES	Human
ACCESSIONS	BC018726
VECTOR	pDONR223
NCBI SYMBOL	CD74
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacaggaggagaagcaggagctgtcgggaagatcagaagccagtcagatgaccagcgcgacattatctccaacaatgagcaactgccatgctggccggcgccctggggcccccggagagcaagtgcagccgggagccctgtacacaggctttccatcctgggtgactctg
ctctcgtggccaggccaccaccgcctacttctgtaccagcagcaggggcgtggacaaactgacagtcacctccagaacctgcagctggagaacctgcgcatgaagcttcccaagcctcccaagcctgtgagcaagatgcgcatggccaccgcgtgatgcaggcgtgc
ccatgggagccctggccaggggccatgcagaatgccaccaagtatggcaaatgacagaggacatgtgatgcacgtgctccagaatgctgacccctgaagggtgtaccggccactgaaggggagcttccggagaacctgagacaccttaagaacacccatggagaccatagact
ggaaggctttgagagctggatgcaccattgctcctgtttgaaatgagcaggcactccttgagcaaaagccactgacgtccaccgaaagagtcactggaactggaggaccgtcttctggcgtgggtgtgaccaagcaggatctggccagctcccatgtga