

BRI3-brain protein I3 Gene

Reference: PTXBC018737

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	brain protein I3; pRGR2		
For research use only.			

GENE DESCRIPTION	brain protein I3
GENE ID	25798
ORIGIN SPECIES	Human
ACCESSIONS	BC018737
VECTOR	pDONR223
NCBI SYMBOL	BRI3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaccacaagccgtctgctgcaggagcggcgccgcctacaacctggaggccggcagggcgactacgcgtcgcccgccgacgggtacgggccatccccgcgcgccccgcgcgcgcctacccctacctcgtcacaggataccaccactatccagggtctacaac
atccacagccggaccgtcaccgctatctgccaactctatcgtggtcgttagggaggtgtctctgctgcagggttgggggtgctggaggactgcttcacctctctgggcatctctctggccatcatctcttcccccttgggttcatttgcgtgtttgccttgaggagcgacgatgccccaaactgtg
gagccaccttgcctaa