

LRRC29-leucine rich repeat containing 29 Gene

Reference: PTXBC018785

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FBL9; FBXL9; leucine-rich repeat-containing protein 29; F-box and leucine-rich repeat protein 9; F-box protein FBL9; F-box/LRR-repeat protein 9; leucine rich repeat containing 29		

For research use only.

GENE DESCRIPTION	leucine rich repeat containing 29
GENE ID	26231
ORIGIN SPECIES	Human
ACCESSIONS	BC018785
VECTOR	pDONR223
NCBI SYMBOL	LRRC29
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtacagctctgggtggcctgcaggagctgcagagcctgcagatggccgagggcgggaactggccagggccctgggctgtatgcacgggctccatccagctggcctccctcagctggccactgctcttcattgaagtcacgccagagctggagcatcaggcctcaggtacca
aggacgctgtccagagccacagggccctccctgctcagctgcgggcccctgcaggagttggactcacagcctgcagcaagctgactgatgccagtttagccaagggtgctccagtttctccagctgaggcagctgtcccttagcctgtgtccagaactcacagacaacgcttggttg
ctgtggccaggggctgtcctagcctggagcactggcgtgagtcactgcagccgactcagtgacaagggtgggcccagggcagccagctcctggccaaggctgcagcatcacaactgtccagctgcagtcagctcatagagcagacactggtatgtatggcaggcgtgcaggca
gctccgggtgttgatgtgccacgtgccctggcatcaacatggccgctcagacgcttccaagccagctgcccaggtgtctgtgtccagtcctcgttgagggggctgacctgacctaacacttga