

WDR82-WD repeat domain 82 Gene

Reference: PTXBC018941

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MST107; MSTP107; PRO2730; PRO34047; SWD2; TMEM113; WDR82A; WD repeat-containing protein 82; transmembrane protein 113; WD repeat domain 82		

For research use only.

GENE DESCRIPTION	WD repeat domain 82
GENE ID	80335
ORIGIN SPECIES	Human
ACCESSIONS	BC018941
VECTOR	pDONR223
NCBI SYMBOL	WDR82
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacctccctgccaccacaaattgtgtacatagcttcagatgataccacccttccccagctccaaccaagagctggttctaggcctgtgttatatgtcatatttagcggttttatatagaccttgattctgtgtttgattttagcacagtgtatgcaccttcatttaatacatctgtgtgcata
cagatacgcatatattgtgtgctgtatgcataatctctcatctgtagttccaagagttcagctgaagcagatggagtcctgcagcccaggagacacctgcacccctgctaatagtgtttgccacaagtatttagtgagttctcttattaattttcatttcagaagactgaagcaagctgatagt
gtttgctgtttctttggcagctaagtga