



Reference: PTXBC018950

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSP 75; HSP75; HSP90L; TRAP-1; heat shock protein 75 kDa, mitochondrial; TNFR-associated protein 1; testicular tissue protein Li 209; tumor necrosis factor type 1 receptor-associated protein; TNF receptor associated protein 1		

For research use only.

GENE DESCRIPTION	TNF receptor-associated protein 1
GENE ID	10131
ORIGIN SPECIES	Human
ACCESSIONS	BC018950
VECTOR	pDONR223
NCBI SYMBOL	TRAP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr