

RPL11-ribosomal protein L11 Gene

Reference: PTXBC018970

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DBA7; GIG34; L11; 60S ribosomal protein L11; CLL-associated antigen KW-12; cell growth-inhibiting protein 34; ribosomal protein L11		

For research use only.

GENE DESCRIPTION	ribosomal protein L11
GENE ID	6135
ORIGIN SPECIES	Human
ACCESSIONS	BC018970
VECTOR	pDONR223
NCBI SYMBOL	RPL11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggatcaagtgaaaaggagaaccccatgcgggaacttcgcatccgcaaaactctgtctcaacatctgtgtgggagagtgagacagactgacgcgagcagccaagggtgtggagcagctcacaggcagacccctgtgtttccaaagctagatacactgtcagatccttggc
atccggagaaatgaaaagattgtctccactgcacagttcgaggggccaaggcagaagaatacttgagaagggtctaaagggtcgaggagatgagtttaagaaaaaacaactctcagatactggaaacttgggtttggatccaggaaacacatcgtctgggtatcaaatgacccaag
cattggtatctacggcctggacttctatgtgtgctggttaggccagggttcagatcgcagacaagaagcgcaggacaggtgcatggggccaaacacagaatcagcaaaaggaggccatgcgctggttcagcagaagtatgatgggatcatcctctggcaataa