

ZNF148-zinc finger protein 148 Gene

Reference: PTXBC018971

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BERF-1; BFCOL1; HT-BETA; ZBP-89; ZFP148; pHZ-52; zinc finger protein 148; CACCC box-binding protein; CLL-associated antigen KW-10; transcription factor ZBP-89; zinc finger DNA-binding protein 89		

For research use only.

GENE DESCRIPTION	zinc finger protein 148
GENE ID	7707
ORIGIN SPECIES	Human
ACCESSIONS	BC018971
VECTOR	pDONR223
NCBI SYMBOL	ZNF148
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgtttacacagagtagctgcttcagagacatgagaagattcactactgtgtgaaaaaccatttcgctgtgatgaatgtggtgatgagacaaccatgcttctgccacatcaacacagagactttctggatcaagtgaacttcagaagaaagctgaggccagcctgtccaccaagcttacca
aafgagctcctttgaacagcccttcgctgctccatcatggatcaagagctggaatagctactcaatttagcactgccaatggacaggtgaaccttcggggaccagggaagaatgctgaatttcagaatttccttggatgaatgataatagagctgggatgacatcttcacctgatgc
cacaactggccagacttttgctaa